

CEKA ATTACHMENTS PRECI-LINE

PERMA-RET



PERMA-RET

MECHANICAL RETENTION FOR ACRYLIC AND COMPOSITE TEETH



1



2




3



4

A recess is drilled with the twist drill.
The precious metal wire is threaded into the
recess and cut. Bend for additional retention.



1601

PERMA-RET KIT:

- 1 wire holder
- 100 cm of PERMA-RET wire
- 3 twist drills

1602

REFILL:

- 100 cm of PERMA-RET wire

WWW.CEKA-PRECILINE.COM

SWISS MADE



CEKA REVAX
CEKA SOL
PRECI-HORIX
PRECI-VERTIX
PRECI-VERTIX AT
PRECI-BAR
PRECI-BALL
PRECI-CLIP
PRECI-SAGIX
PRECI-CLIX AXIAL
PRECI-CLIX EXTRACORONAL
PRECI-CLIX RADICULAR
PRECI-CLIX RADICULAR RC
PRECI-TUBIX/MORTIX
PRECI-POST
PRECI-PROFILE
PERMA-RET
3C-BOND
EXPANDO
CEKA-MAX

PERMA-RET

Prospectus no 260.E-ed.C
Copyright© ALPHADENT NV 8/2010
BTW BE 0414.789.321 RPR Kortrijk

ALPHADENT NV
Textielstraat 24
8790 Waregem, Belgium
tel +32 56 629 100
fax +32 56 629 101
e-mail: info@ceka-preciline.com